

## Recombinant SARS-CoV-2 S1, 5 mg

---

<b>Catalog number</b>	710092
<b>Description</b>	Recombinant human coronavirus SARS-CoV-2 Spike S1 protein contains amino acids 14-681 of the S protein, and two extra amino acids (AS) at the N-terminus and a His-6 tag at the C-terminus. Predicted molecular weight: 76.4 kDa.
<b>Amino acid sequence</b>	-
<b>Product host</b>	CHO cells
<b>Product formulation</b>	Liquid
<b>Product buffer solution</b>	PBS, pH 7.4
<b>Reconstitution</b>	-
<b>Shelf life and storage</b>	Unspecified, storage at -70 °C
<b>Analyte description</b>	Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) is an enveloped, positive stranded RNA virus. SARS-CoV-2 causes respiratory illness called COVID-19.
<b>Product concentration</b>	0.5–1 mg/ml
<b>Purity</b>	SDS-PAGE > 90 %
<b>Antibodies tested</b>	-

**Legal disclaimer**