

Tricyclic Antidepressant (TCA) Antibody

Product overview

Catalog number	100896
Specificity	Antibody recognizes Tricyclic Antidepressants
Description	Protein A purified monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components
Product buffer solution	0.05M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl, 0.05% NaN ₃ as a preservative
Shelf life and storage	18 months from manufacturing at 2–8 °C
Subclass	IgG ₁
Analyte description	Tricyclic antidepressants are a class of medications used primarily for the treatment of depression and other mood disorders.

Parameters tested on each lot

Product appearance	Liquid, may turn slightly opaque during storage
Product concentration	≥ 1.5 mg/ml
Immunoreactivity	80–120 % compared to the reference sample in an FIA test
IEF Profile	6.1–7.1
Purity	≥ 95 %

Kinetic parameters

Association rate constant	N/D
Dissociation rate constant	N/D
Affinity constant	N/D
Determination method	-
Determination antigen	-



Legal disclaimer

Cross-reactivities	N/D				
Epitope	N/D				
Pair recommendations	T001A or CB500 (BSA conjugate as capture) Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.				
Platforms tested	ELISA, LF				
Antigens tested	Please see Pair recommendations above				
Product stability	<table><thead><tr><th>TEMPERATURE, TIME</th><th>RESULT</th></tr></thead><tbody><tr><td>N/D</td><td>-</td></tr></tbody></table>	TEMPERATURE, TIME	RESULT	N/D	-
TEMPERATURE, TIME	RESULT				
N/D	-				
Miscellaneous	This product was previously available with Catalog codes T002M and CB300. It has been produced using the same original mouse monoclonal hybridoma cell line.				
References	-				



Legal disclaimer