

Zolpidem-BSA Antigen

Catalog number	170041
Description	Synthetic Zolpidem, BSA conjugated
Amino acid sequence	-
Product host	-
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	Shelf life unspecified. Short term storage (≤ 2 weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Analyte description	Zolpidem is a central nervous system depressant drug causing drowsiness, sedation and decreased psychomotor function. Zolpidem is used to treat insomnia and may result in next-day impairment.
Product concentration	Lot-specific
Matched pairs	P03-99-11R-IF, P03-99-12R-IF, P03-99-13R-IF, JP046, JP047, JP048, HM311 or HM312 (detection)
Miscellaneous	This product was previously available with a Catalog code P83-99-11A. It has been produced using the same conjugation protocol.

**Legal disclaimer**