

## Zolpidem-BSA Antigen

---

<b>Catalog number</b>	170041
<b>Description</b>	Synthetic Zolpidem, BSA conjugated
<b>Amino acid sequence</b>	-
<b>Product host</b>	-
<b>Product formulation</b>	Liquid; Antigen-BSA conjugate
<b>Product buffer solution</b>	Distilled water
<b>Reconstitution</b>	-
<b>Shelf life and storage</b>	60 months from manufacturing at -20 °C. Short term storage ( $\leq 2$ weeks): 2–8 °C. Avoid repeated freezing and thawing.
<b>Analyte description</b>	Zolpidem is a central nervous system depressant drug causing drowsiness, sedation and decreased psychomotor function. Zolpidem is used to treat insomnia and may result in next-day impairment.
<b>Product concentration</b>	Lot-specific
<b>Matched pairs</b>	JP046, JP047, JP048, HM311 or HM312 (detection)
<b>Miscellaneous</b>	This product was previously available with a Catalog code P83-99-11A. It has been produced using the same conjugation protocol.

**Legal disclaimer**