

PRODUCT SPECIFICATIONS

Name	Anti-h FABP3 2302 SPTN-5
Specificity	Antibody recognizes human fatty acid binding protein 3 (FABP3)
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components
Product code	100292
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Shelf life and storage	18 months from manufacturing at 2-8 °C
Analyte description	The fatty-acid binding proteins (FABPs) are a family of carrier proteins for fatty acids and other lipophilic substances. Heart-type fatty acid-binding protein (H-FABP= FABP3) is a low molecular weight cytoplasmic protein present abundantly in the myocardium. When the myocardium is injured, as in the case of myocardial infarction, low molecular weight cytoplasmic proteins including H-FABP are released into the circulation and H-FABP can be used as a biochemical diagnostic marker in the early phase of acute myocardial infarction (AMI).

PARAMETERS TESTED ON EACH LOT

Product appearance	Liquid, may turn slightly opaque during storage
Product concentration	5.0 mg/ml (+/- 10 %)
Immunoreactivity	80-120 % compared to the reference sample in an FIA test
IEF Profile	6.9 - 7.6
Purity	≥ 95 %

PARAMETERS DETERMINED DURING PRODUCT DEVELOPMENT

Subclass	IgG ₁
Association rate constant	2.9×10^5 1/Ms
Dissociation rate constant	4.5×10^{-3} 1/s
Affinity constant	$K_A = 6.5 \times 10^7$ 1/M; $K_D = 1.6 \times 10^{-8}$ M (= 16 nM)
Determination method	SPR analysis (ProteOn XPR36)
Determination antigen	FABP3, Lee Biosolutions (Cat 276-10, Lot L84659-A)

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

Oy Medix Biochemica Ab

Headquarters and manufacturing
Klovinpellontie 3
FI-02180 Espoo, Finland

medix@medixbiochemica.com
www.medixbiochemica.com
VAT reg. no. FI4631532

Cross-reactivities No cross-reactivity against FABP1 or FABP4, others not tested.

Epitope Not Determined (N/D)

Pair recommendations	CAPTURE ANTIBODY	DETECTION ANTIBODY
	2302, 2303	2304
	2304	2302, 2303

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Product stability	TEMPERATURE, TIME	RESULT
	-70 °C, 21 days	OK
	-20 °C, 21 days	OK
	+4 °C, 21 days	OK
	+25 °C, 21 days	OK
	+35 °C, 21 days	OK
	+45 °C, 7 days	OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2-8 °C in the product buffer.

Miscellaneous -

References -

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

Oy Medix Biochemica Ab

Headquarters and manufacturing
Klovinpellontie 3
FI-02180 Espoo, Finland

medix@medixbiochemica.com
www.medixbiochemica.com
VAT reg. no. FI4631532