

PRODUCT SPECIFICATIONS

Name	Anti-h Transferrin 3106 SPTN-5
Specificity	Antibody recognizes human transferrin
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components
Product code	100474
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Shelf life and storage	Unspecified, storage at 2-8 °C
Analyte description	Transferrins are glycoproteins which regulate the amount of free iron in circulation. Each transferrin molecule binds two Fe ⁺³ with high affinity. When no iron is bound the molecule is called Apotransferrin. Transferrin concentration elevates in anemia, but also during pregnancy and when using birth control pills. Chronic infections, malignant growth, alcohol induced liver damage, and other conditions with protein loss may lower transferrin concentration.

PARAMETERS TESTED ON EACH LOT

Product appearance	Liquid, may turn slightly opaque during storage
Product concentration	5.0 mg/ml (+/- 10 %)
Immunoreactivity	80-120 % compared to the reference sample in an FIA test
IEF Profile	7.4 - 8.7
Purity	≥ 95 %

PARAMETERS DETERMINED DURING PRODUCT DEVELOPMENT

Subclass	IgG _{2a}
Association rate constant	2.2 x 10 ⁶ 1/Ms
Dissociation rate constant	1.1 x 10 ⁻⁴ 1/s
Affinity constant	K _A = 2.1 x 10 ¹⁰ 1/M; K _D = 4.8 x 10 ⁻¹¹ M (=0.048 nM)
Determination method	SPR analysis (ProteOn XPR36)
Determination antigen	Transferrin native human, AbD Serotec (Cat PHP-238, Lot 171210)

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

Oy Medix Biochemica Ab

Headquarters and manufacturing
Klovinpellontie 3
FI-02180 Espoo, Finland

medix@medixbiochemica.com
www.medixbiochemica.com
VAT reg. no. FI4631532

Cross-reactivities Not determined (N/D)

Epitope N/D

Pair recommendations	CAPTURE ANTIBODY	DETECTION ANTIBODY
	3104, 3105	3106
	3106	3104, 3105

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Product stability	TEMPERATURE, TIME	RESULT
	-70 °C, 21 days	OK
	-20 °C, 21 days	OK
	+4 °C, 21 days	OK
	+35 °C, 7 days	OK
	+35 °C, 21 days	not recommended
	+45 °C, 3 days	not recommended

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2-8 °C in the product buffer.

Miscellaneous -

References -

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

Oy Medix Biochemica Ab

Headquarters and manufacturing
Klovinpellontie 3
FI-02180 Espoo, Finland

medix@medixbiochemica.com
www.medixbiochemica.com
VAT reg. no. FI4631532