

PRODUCT SPECIFICATIONS

Name	Recombinant HCoV NL63 N, 100 µg
Description	Recombinant human coronavirus NL63 nucleoprotein with an N-terminal histidine tag. Predicted molecular weight: 42 kDa.
Production host	<i>Escherichia coli (E.coli)</i>
Product code	610041
Product formulation	Lyophilized from 0.2 µm filtered solution
Product buffer	50 mM Tris-HCl, pH 7.5, 500 mM NaCl, 0.1 % CHAPS; 6 % sucrose as stabilizer
Reconstitution	Reconstitute lyophilized protein with 200 µl of deionized water
Shelf life and storage	Shelf life unspecified, storage at 2-8 °C
Analyte description	Coronavirus is believed to cause a significant percentage of all common colds in human adults. Coronaviruses primarily infect the upper respiratory and gastrointestinal tract.
Product concentration	0.5 mg/ml when reconstituted with 200 µl of deionized water
Reactivity with MedixMABs	Anti-hCoV 3501: +

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

Oy Medix Biochemica Ab

Headquarters and manufacturing
Klovinpellontie 3
FI-02180 Espoo, Finland

medix@medixbiochemica.com
www.medixbiochemica.com
VAT reg. no. FI4631532