


PRODUCT SPECIFICATIONS

Name	Recombinant HCoV OC43 N, 100 µg
Description	Recombinant human coronavirus OC43 nucleoprotein with an N-terminal histidine tag. Predicted molecular weight: 49 kDa.
Production host	<i>Escherichia coli</i> (<i>E. coli</i>)
Product code	610040
Product formulation	Lyophilized from 0.2 µm filtered solution
Product buffer	50 mM Tris-HCl, pH 7.5, 500 mM NaCl, 0.1 % CHAPS; 6 % sucrose as a stabilizer
Reconstitution	Reconstitute lyophilized protein with 200 µl of deionized water
Shelf life and storage	Shelf life unspecified, storage at 2-8 °C
Analyte description	Human coronavirus OC43 (hCoV-OC43) is an enveloped, positive-stranded RNA virus. HCoV-OC43 causes a significant fraction of upper respiratory tract infections, especially common colds.
Product concentration	0.5 mg/ml when reconstituted with 200 µl of deionized water
Purity	SDS-PAGE 
Reactivity with MedixMABs	Anti-hCoV 3503: + Anti-hCoV 3504: + Anti-hCoV 3505: + Anti-hCoV 3506: +

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

Oy Medix Biochemica Ab

Headquarters and manufacturing
 Klovinpellontie 3
 FI-02180 Espoo, Finland

medix@medixbiochemica.com
 www.medixbiochemica.com
 VAT reg. no. FI4631532