


PRODUCT SPECIFICATIONS

Name	Recombinant cTnT antigen, 100 µg
Description	Recombinant human cardiac troponin T antigen (isoform 6) with a C-terminal histidine tag. Predicted molecular weight: 36 kDa.
Amino acid sequence	MSDIEEVVEEYEEEEQEEAAVEEQEEAAEEDAEEAEAEETEETRAEEDDEEEEAKEAEDGPMEEESKPKPRSF MPNLVPPKIPDGERVDFDDIHRKRMEKDLNELQALIEAHFENRKKKEEELVSLKDRIERRRAERAEQQ RIRNEREKERQNRLAEERARREEEENRRKAEDEARKKKALSNNMMHFGGYIQKQAQTERKSGKRQTE REKKKKILAERRKVLAIIDLHNLNEDQLREKAKELWQSIYNLEAEKFDLQEKFKQKYEINVLNRINDNQ KYSKTRGKAKVTGRWKLEHHHHHH
Production host	<i>Escherichia coli</i> (<i>E. coli</i>)
Product code	610101
Product formulation	Lyophilized
Product buffer	50 mM Tris-HCl, pH 7.5; 150 mM NaCl, 4 M urea; containing 3 % sucrose, 2 % D-mannitol, and 0.01 % Tween 20 as stabilizers
Reconstitution	Reconstitute the lyophilized protein with 200 µl of deionized water
Shelf life and storage	Unspecified, storage at 2-8 °C
Analyte description	Troponins form a complex of three regulatory proteins (troponin C, I and T) that are integral to muscle contraction in skeletal and cardiac muscles. Serum cardiac troponin tests can be used to help diagnose several different heart disorders, especially myocardial infarction.
Product concentration	0.5 mg/ml when reconstituted with 200 µl of deionized water
Purity	SDS-PAGE 
Reactivity with MedixMABs	Not applicable (N/A)

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

Oy Medix Biochemica Ab

Headquarters and manufacturing
Klovinpellontie 3
FI-02180 Espoo, Finland

medix@medixbiochemica.com
www.medixbiochemica.com
VAT reg. no. FI4631532