


PRODUCT SPECIFICATIONS

Name	Recombinant Influenza A NP, 100 µg
Description	Recombinant influenza A nucleoprotein antigen (strain A/Hong Kong/1/1968 H3N2) with a histidine tag. Predicted molecular weight: 57 kDa
Production host	<i>Escherichia coli (E. coli)</i>
Product code	610050
Product formulation	Liquid
Product buffer solution	50 mM Tris-HCl, 250 mM NaCl, 10 % glycerol, pH 8.0
Shelf life and storage	Unspecified, storage at -20 °C
Analyte description	Influenza A is a genus of the Orthomyxoviridae family of viruses. It includes one viral species which causes influenza in birds and some mammals. Influenza A subtypes are labeled according to their hemagglutinin and neuraminidase protein variants. Influenza nucleoprotein (NP) binds to viral RNA, and these ribonucleoprotein complexes play a major role in viral transcription and replication. NP is highly conserved among influenza A virus strains, making NP a very suitable control antigen to be used in influenza A immunological tests.
Product concentration	1.0 mg/ml (+/- 10 %)
Purity	SDS-PAGE 
Reactivity with MedixMABs	Anti-Influenza A 7304: + Anti-Influenza A 7305: + Anti-Influenza A 7307: + Anti-Influenza A 7314: +

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks or registered trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

Oy Medix Biochemica Ab

Headquarters and manufacturing
 Klovinpellontie 3
 FI-02180 Espoo, Finland

medix@medixbiochemica.com
 www.medixbiochemica.com
 VAT reg. no. FI4631532