

PRODUCT SPECIFICATIONS

Name	Vitamin D antigen, BSA conjugated, 100 µg
Description	Synthetic Vitamin D antigen, BSA conjugated. The conjugation position is C22.
Amino acid sequence	-
Production host	-
Product code	710060
Product formulation	Liquid, turbid
Product buffer solution	0.9 % NaCl
Reconstitution	-
Shelf life and storage	60 months, short time storage at -20 °C, long time storage at -80 °C
Analyte description	Vitamin D is a fat-soluble vitamin that is made by the body after exposure to ultraviolet rays. Vitamin D increases intestinal adsorption of calcium, magnesium and phosphate and is needed for normal bone growth. In humans, Vitamin D occurs in several forms including vitamin D2 and D3.
Product concentration	> 0.5 mg/ml
Purity	Not Determined (N/D)
Reactivity with MedixMABs	3850, 3851
Other	Reactivity with other supplier's antibodies is not probable due to the unique conjugation position.

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks or registered trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

Oy Medix Biochemica Ab

Headquarters and manufacturing
Klovinpellontie 3
FI-02180 Espoo, Finland

medix@medixbiochemica.com
www.medixbiochemica.com
VAT reg. no. FI4631532