

## PRODUCT SPECIFICATIONS

<b>Name</b>	Anti-h GH 5802 SP-1
<b>Specificity</b>	Antibody recognizes human growth hormone (GH)
<b>Description</b>	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components
<b>Product code</b>	100042
<b>Product buffer solution</b>	0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
<b>Shelf life and storage</b>	36 months from manufacturing at 2-8 °C
<b>Analyte description</b>	Growth hormone (GH) is a single-chain polypeptide hormone that is synthesized, stored, and secreted by the anterior pituitary gland. It stimulates growth and cell reproduction and regeneration in humans and other animals. Growth hormone is used clinically to treat children's growth disorders and adult growth hormone deficiency.

## PARAMETERS TESTED ON EACH LOT

<b>Product appearance</b>	Liquid, may turn slightly opaque during storage
<b>Product concentration</b>	1.0 mg/ml (+/- 10 %)
<b>Immunoreactivity</b>	80-120 % compared to the reference sample in an FIA test
<b>IEF Profile</b>	6.0 – 7.1
<b>Purity</b>	≥ 95 %

## PARAMETERS DETERMINED DURING PRODUCT DEVELOPMENT

<b>Subclass</b>	IgG <sub>1</sub>
<b>Association rate constant</b>	Not Determined (N/D)
<b>Dissociation rate constant</b>	N/D
<b>Affinity constant</b>	1 x 10 <sup>10</sup> l/mol
<b>Determination method</b>	Radioimmunoassay (RIA)
<b>Determination antigen</b>	GH, Boehringer (Cat 252964, Lot 1222303)

### Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB is a registered trademark of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

### Oy Medix Biochemica Ab

**Headquarters and manufacturing**  
Klovinpellontie 3  
FI-02180 Espoo, Finland

medix@medixbiochemica.com  
www.medixbiochemica.com  
VAT reg. no. FI4631532

**Cross-reactivities**

PRL 1.0 % (International Laboratory Service, Lot 11385)  
 PL <0.2 % (UCB-Bioproducts, Cat i020, Lot 007)

**Epitope**

N/D

**Pair recommendations**

**CAPTURE ANTIBODY**  
 5802

**DETECTION ANTIBODY**  
 5801

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

**Product stability**

TEMPERATURE, TIME	RESULT
-70 °C, 21 days	N/D
-20 °C, 21 days	N/D
+4 °C, 21 days	N/D
+35 °C, 7 days	N/D
+35 °C, 21 days	N/D
+45 °C, 3 days	N/D
+45 °C, 7 days	N/D

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2-8 °C in the product buffer.

**Miscellaneous**

The antibody recognizes growth hormone isoform 22K, and does not recognize isoform 20K. Others not tested.

**References**

Erb, T., Karolyi, G., Pampalone, A., Eberle, A.N., Baumann, J.B., Bürgisser, E., and Girard, J., (1988) Human growth hormone in urine: development of an ultrasensitive radiometric assay. *Experientia*, 44:155-157