



We are Actim

Leading the way to right
treatment.



Actim is developing three different diagnostic tests to aid in COVID-19 pandemic

Available for evaluation

Serological antibody detection test **ELISA SARS-CoV-2 IgG**

- An immunoenzymometric assay for quantitative detection of specific IgG to SARS-CoV-2
- Sample type: serum from venous blood
- The test is intended for professional use to help diagnose recent COVID-19 infection
- Plates and reagents are developed and produced in Finland

Upcoming

Serological antibody detection test **Actim SARS-CoV-2 IgG/IgM**

- A qualitative immunochromatographic dipstick test for detection of specific IgG and IgM to SARS-CoV-2
- Sample type: small amount of whole blood
- The test is intended for professional use to help diagnose recent COVID-19 infection
- Tests are developed and produced in Finland

Upcoming

SARS-CoV-2 virus antigen detection rapid test for acute use **Actim SARS-CoV-2**

- A qualitative immunochromatographic dipstick test for detection of SARS-CoV-2 virus antigens
- Sample type: respiratory sample
- The test is intended for professional use to help diagnose acute COVID-19 infection
- Tests are developed and produced in Finland