

Product specifications

Name	Vitamin D antigen, Biotin conjugated, 1 mg
Description	Synthetic Vitamin D antigen, Biotin conjugated.
Amino acid sequence	-
Product host	-
Product code	710049
Product formulation	Lyophilized antigen, on the bottom of the vial. Not visually detectable.
Product buffer solution	-
Reconstitution	Reconstitute the lyophilized antigen with 1000 µl of EtOH
Shelf life and storage	Lyophilized product 60 months at -20 °C. After reconstitution 6 months at -20 °C. Long term storage at -80 °C.
Analyte description	Vitamin D is a fat-soluble vitamin that is made by the body after exposure to ultraviolet rays. Vitamin D increases intestinal adsorption of calcium, magnesium and phosphate and is needed for normal bone growth. In humans, Vitamin D occurs in several forms including vitamin D2 and D3.
Product concentration	1 mg/ml when reconstituted with 1000 µl of EtOH
Purity	Not Determined (N/D)
Reactivity with MedixMABs	Anti-Vitamin 25 OH D2/D3 3850: + Anti-Vitamin 25 OH D2/D3 3851: +
Other	The synthetic antigen is designed to perform as a matching pair with the antibody 3850 or 3851. Its reactivity with antibodies from other suppliers cannot be guaranteed due to the unique conjugation position.

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

