Oy Medix Biochemica Ab
Headquarters
Klovinpellontie 3
FI-02180 Espoo
Finland
Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.
10850 Metro Court
Maryland Heights, MO 63043
United States of America
Email: sales@medixbiochemica.com

Product Datasheet

Vaginal Swabs, Single Donor, Pre-Coitus
991-25-PREC

Product Name	Vaginal Swabs, Single Donor, Pre-Coitus
Catalog Number	991-25-PREC
Description	Each sample is equivalent to 2 swabs per donor. Swabs are taken a minimum of 6 days following last intercourse. Inquire for larger donor counts, swab counts, or custom criteria.
Tested Applications	QC Materials
Alternative Names	Vaginal, vagina, reproductive, reproduction, genital, vaginal culture, gynecological, intercourse, sex, sexual
Brand	Lee Biosolutions
Form/Appearance	Swab
Storage	-20ºC
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None added, inquire if needed